

Sequencing Facility

CCR-Sequencing Facility Illumina Sequencing Report

Project Information

Principal Investigator: Jing Wu

PI Laboratory Contact: Madison Butler

Bioinformatics Contact: CCBR, Parthav Jailwala

Project Title: JingWu_CS026381_44RNA_021320

NAS Order ID: CS026381

Samples Total in project: 44

Samples in This Report: 44

Completion of NAS: Yes

Report Date: 03/14/20

Sequencing Details

Flowcell ID:	HJJKWDRXX	Sequence Control:	PhiX
Instrument:	NovaSeq	Control Result:	Pass
Sequencing Type:	mRNA-Seq	Library Protocol:	TruSeq Stranded mRNA Library Prep
Read Length:	151 (2x151 cycles)	Sequencing Chemistry:	NovaSeq
Multiplexing:	44 per lane	Reference Genome:	Human_hg38
Strand Specificity:	Stranded	Annotation:	Ensembl96_30 GTF

Run Comments

44 mRNA-Seq samples were pooled and sequenced on NovaSeq_SP using Illumina TruSeq Stranded mRNA Library Prep and paired-end sequencing. The samples have 33 to 50 million pass filter reads with more than 94% of bases above the quality score of Q30. Reads of the samples were trimmed for adapters and low-quality bases using Cutadapt before alignment with the reference genome (Human_hg38) and the annotated transcripts using STAR. The average mapping rate of all samples is 95%. Unique alignment is above 90%. There are 3.35 to 5.93% unmapped reads. The mapping statistics are calculated using Picard software. The samples have 0.00% ribosomal bases. Percent coding bases are between 54-60%. Percent UTR bases are 32-35%, and mRNA bases are between 89-92% for all the samples. Library complexity is measured in terms of unique fragments in the mapped reads using Picard's MarkDuplicate utility. The samples have 72-77% non-duplicate reads. In addition, the gene expression quantification analysis was performed for all samples using STAR/RSEM tools. Both the normalized count and the raw count are provided as part of the data delivery.

Note: Residual samples will be retained up to **90 days** of the delivery of this report. To avoid shipping charges, please contact SFILLUMINALAB@mail.nih.gov to arrange pickup samples prior to this time.

Note: Sequencing data will be available to download for **two weeks** following delivery of this report. Please download the data files as soon as possible.


For questions on any aspect of this report please contact CCRSF_IFX@nih.gov.



Leidos Biomedical Research, Inc.

<https://ostr.cancer.gov/resources/fnl-cores/sequencing-facility>

Analysis Workflow



Software and Parameters


Analysis Step	Software	Software Parameters / Notes
Basecalling	RTA v3.4.4	Illumina instrument run time analysis software
Demultiplexing	Bcl2fastq v2.17	--no-lane-splitting -i RunFolder/Data/Intensities/BaseCalls -R RunFolder -barcode-mismatches 1 --ignore-missing-bcls --ignore-missing-filter --ignore-missing-positions --ignore-missing-controls --sample-sheet SampleSheet.csv -o Unaligned
Filtering (Adaptor and quality)	Cutadapt 1.18	-j 8 -b file:adapters.fa -B file:adapters.fa --nextseq-trim=2 --trim-n -n 5 -O 5 -q 10,10 -m 35:35 -o trimmed_R1.fq -p trimmed_R2.fq input_R1.fq input_R2.fq
Alignment	STAR 2.7.0f	<p>1-pass: --genomeDir \$star_genome --outSAMunmapped Within --outFilterType BySJout --outFilterMultimapNmax 20 --outFilterMismatchNmax 999 --outFilterMismatchNoverLmax 0.04 --alignIntronMin 20 --alignIntronMax 1000000 --alignMatesGapMax 1000000 --alignSJoverhangMin 8 --alignSJDBoverhangMin 1 --sjdbScore 1 --readFilesCommand zcat --readFilesIn \$trimmed_R1.fastq.gz \$trimmed_R2.fastq.gz --runThreadN numThreads --outFilterMatchNminOverLread 0.66 --outSAMtype BAM Unsorted --quantMode TranscriptomeSAM --peOverlapNbasesMin 10 --alignEndsProtrude 10 ConcordantPair</p> <p>2-pass: --genomeDir \$star_genome --outSAMunmapped Within --outFilterType BySJout --outFilterMultimapNmax 20 --outFilterMismatchNmax 999 --outFilterMismatchNoverLmax 0.04 --alignIntronMin 20 --alignIntronMax 1000000 --alignMatesGapMax 1000000 --alignSJoverhangMin 8 --limitSjdbInsertNsj 2500000 --sjdbFileChrStartEnd \$input_1-path_sj --alignSJDBoverhangMin 1 --sjdbScore 1 --readFilesCommand zcat --readFilesIn \$trimmed_R1.fastq.gz \$trimmed_R2.fastq.gz --runThreadN \$numthreads --outFilterMatchNminOverLread 0.66 --outSAMtype BAM Unsorted --quantMode TranscriptomeSAM --peOverlapNbasesMin 10 --alignEndsProtrude 10 ConcordantPair</p>

For questions on any aspect of this report please contact CCRSF_IFX@nih.gov.

Sequencing Facility

RNAStatistics	Picard 2.18.26	CollectRnaSeqMetrics.jar REF_FLAT=annotation_refFlat.txt INPUT=sample.bam OUTPUT= RnaSeqMetrics.txt RIBOSOMAL_INTERVALS= ribosome_interval_list.txt STRAND_SPECIFICITY=SECOND_READ_TRANSCRIPTION_STRAND VALIDATION_STRINGENCY=LENIENT
Duplication Statistics	Picard 2.18.26	MarkDuplicates.jar INPUT=sample.bam OUTPUT=sample.MKDUP.bam METRICS_FILE=sample.bam.metric ASSUME_SORTED=true MAX_FILE_HANDLES_FOR_READ_ENDS_MAP=1000 VALIDATION_STRINGENCY=LENIENT
Quantification	RSEM 1.3.1	rsem-calculate-expression -bam --paired-end --estimate-rspd Transcriptome.out.bam \$RSEM_Genome \$Sample_Name

Data Statistics




For questions on any aspect of this report please contact CCRSF_IFX@nih.gov.



Leidos Biomedical Research, Inc.

<https://ostr.cancer.gov/resources/fnl-cores/sequencing-facility>

Sequencing Facility




For questions on any aspect of this report please contact CCRSF_IFX@nih.gov.



Leidos Biomedical Research, Inc.

<https://ostr.cancer.gov/resources/fnl-cores/sequencing-facility>

RNA Statistics



For questions on any aspect of this report please contact CCRSF_IFX@nih.gov.



Leidos Biomedical Research, Inc.

<https://ostr.cancer.gov/resources/fnl-cores/sequencing-facility>

Sequencing Facility

Notes

- **Sample Yield** – The sum of all bases in reads that passed filtering per sample. Indicates the output in million bases (Mb) per lane.
- **% \geq Q30** – The percentage of bases called with an inferred accuracy of 99.9% or above, a measure of basecalling quality.
- **% Total (Primary) Alignment** – The percentage of filtered reads that align to the reference; for mRNA-seq, to the reference genome and the splice junctions. Reads aligning to multiple locations are included in the calculation
- **% Unique Alignment** – The percentage of filtered reads that align uniquely to the reference; for mRNA-Seq, the reference genome and known splice junctions. Reads aligning to multiple locations and abundant sequences are not included in the score.
- **% Non-duplicated Reads** – The percentage of aligned reads with non-redundant start coordinate.
- **% RNA Statistics** – Collect metrics about the alignment of RNA to various functional classes of loci in the genome: coding, intronic, UTR, intergenic, ribosomal. Also determines strand-specificity for strand-specific libraries.

PCT_RIBOSOMAL_BASES: RIBOSOMAL_BASES / PF_ALIGNED_BASES

PCT_CODING_BASES: CODING_BASES / PF_ALIGNED_BASES

PCT_UTR_BASES: UTR_BASES / PF_ALIGNED_BASES

PCT_INTRONIC_BASES: INTRONIC_BASES / PF_ALIGNED_BASES

PCT_INTERGENIC_BASES: INTERGENIC_BASES / PF_ALIGNED_BASES

PCT_MRNA_BASES: PCT_UTR_BASES + PCT_CODING_BASES

For questions on any aspect of this report please contact CCRSF_IFX@nih.gov.



Leidos Biomedical Research, Inc.

<https://ostr.cancer.gov/resources/fnl-cores/sequencing-facility>